



APZEM PREMIUM BACK DRAFT GROSSING STATION (Rbi Series)

APZEM Premium series Grossing Station were designed to provide Sturdy, customizable and affordable workstation for medium to large Labs. The premium series workstation comes with higher suction capacity and additional features that assists in grossing as well as protects you from inhaling harmful formalin and other fumes while grossing.



MODEL NO	SIZE (LENGTH)
RB500i	5 FEET
RB600i	6 FEET
RB700i	7 FEET



STANDARD DESIGN FEATURES

- ▶ All Stainless steel Construction made of 16 Guage SS 304.
- ▶ Inbuilt ventilation system capable of removing all formalin fumes from work top and from Technician's breathing zone.
- ▶ Deep Sink with Cold water Faucet.
- ▶ Back Draft fume Removal system.
- ▶ Removable perforated worktop.
- ▶ Powerful LED Lighting.
- ▶ C-Fold paper towel Holder.
- ▶ Polyethylene Dissecting Board for Dissecting and for Photography.
- ▶ Magnetic Tool holder.
- ▶ Cleanable Pre Filter to remove fine particles.
- ▶ Worktable Height : 800mm.
- ▶ Additional shelf for storing.
- ▶ Hand rinsing assembly.
- ▶ Perimeter Rinse assembly.
- ▶ Ventilated Trash cabinet.
- ▶ High capacity 0.55 HP waste disposal unit.
- ▶ Spot extraction system.
- ▶ Transparent Lexan air deflector panels.
- ▶ 40mm height adjustable leg for accurate leveling of work surface.
- ▶ Built in storage cabinet.
- ▶ Formalin Dispensing system.

OPTIONAL FEATURES AND ACCESSORIES

- ▶ Magnifying Lamp fixture for mounting 3 axis adjustable magnifying Lamp.
- ▶ Digital camera holder.
- ▶ Safety splash shield.
- ▶ CPU Mounting bracket.
- ▶ Mobile storage cabinet with draw arrangement for storing.
- ▶ Adjustable Shelving.
- ▶ Compact storage box with lock and key.
- ▶ Higher grade, SS 316 Construction for better chemical resistance and greater life.
- ▶ Optional Carbon filter / potassium permanganate filter.
- ▶ Additional self contained exhaust blower for Increased suction.
- ▶ Pull out writing pad with additional storage space to store documents.

** Chose any of 48 accessories listed at end..