




APZEM ELEVATING DOWNDRAFT GROSSING STATION

(DHE SERIES)

APZEM Elevating series downdraft Grossing Station were designed to provide Sturdy, customizable and affordable workstation for medium to large Labs. The Elite Elevating series workstation comes with higher suction capacity, ergonomically designed and additional height adjustable features that assists in high productivity grossing as well as protects you from inhaling harmful Formalin fumes.



MODEL NO	SIZE (LENGTH)
DHE6000	6 FEET

Specification

Model	DH6000E
MOC	SS304 (Standard) / SS316 (Optional)
Ventilation	Downdraft
Dimension	1800mm Length x 905mm Width x 1885mm Height
Work Surface	900mm height from the ground Height can be increased upto 1200mm
Sink	790mm x 660mm x 310mm deep
Filters	Water washable polyester prefilters along with optional custom designed filter
Electrical	220-240 VAC 50/60 Hz, 1 Phase, 110 V supplied on request
Plumbing	1/2" industrial grade flexible hose supplied for hot and cold water connection 1 1/2" drain hose (waste disposal)
Waste Disposal	0.55 HP High capacity waste disposal unit.
Exhaust air capacity	1800 CMH
Perforated working area	800mm length x 750mm width
exhaust duct diameter	200mm diameter

Standard Features



Anti Corrosive Construction



Perforated Worktop



Inbuilt Ventilation



Spacious Sink



Downdraft Extraction



Powerful Lights



Boundry Rinse



Towel Holders




Laser Etched Cutting Board




Instrument Holder (Magnetic)



Path Imaging System



Ergonomic Design



Height Adjustable



Spacious Shelves



Hand Rinse



Bio Decomposer



Transparent Panels



Work floor Levelling



Formalin Dispenser



Formalin Monitoring System

Optional Features



Magnifying Lamp



MOC: SS316



Optional Carbon Filter



Organizer Bin



Digital Weighing Scale



Ruler (Inch scale)