




# APZEM MINI BACKDRAFT GROSSING STATION

(FB SERIES)

APZEM Premium series Backdraft Grossing Station were designed to provide Sturdy, customizable and affordable workstation for medium to large Labs. The premium series workstation comes with higher suction capacity and additional features that assists in grossing as well as protects you from inhaling harmful Formalin and other harmful fumes while grossing.



| MODEL NO | SIZE (LENGTH) |
|----------|---------------|
| FB400    | 4 FEET        |
| FB500    | 5 FEET        |

## Specification

|                         |   |
|-------------------------|---|
| Model                   | FB400   |
| MOC                     | SS304 (Standard) / SS316 (Optional)   |
| Ventilation             | Backdraft   |
| Dimension               | 1275mm Length x 800mm Width x 1960mm Height   |
| Work Surface            | 900mm height  |
| Sink                    | 380mm length x 350mm width x 180mm deep sink  |
| Filters                 | Water washable polyester prefilters   |
| Electrical              | 220 - 240 VAC 50/60 Hz, 1 Phase, 110 V supplied on request  |
| Plumbing                | 1/2" industrial grade flexible hose supplied for water connection<br>1 1/2" drain hose (waste disposal) |
| Waste Disposal          | 0.55 HP High capacity waste disposal unit.  |
| Exhaust air capacity    | 930 CMH   |
| Perforated working area | 880mm length x 495mm width  |
| Exhaust duct diameter   | 150MM diameter  |

## Standard Features



Anti Corrosive Construction



Perforated Worktop



Inbuilt Ventilation



Spacious Sink



Backdraft Extraction



Powerful Lights



Boundry Rinse



Towel Holders




Laser Etched Cutting Board



Instrument Holder (Magnetic)



Custom Filters



Ergonomic Design



Spacious Shelves



Hand Rinse



Bio Decomposer



Transparent Panels



Work floor Levelling



Camera Holder

## Optional Features



Magnifying Lamp



Safety Shield



MOC: SS316



Optional Carbon Filter



Formalin Monitoring System



Adjustable Shelves



Organizer Bin