




APZEM PROFESSIONAL BACKDRAFT GROSSING STATION

(RB SERIES)

APZEM Premium series Backdraft Grossing Station were designed to provide Sturdy, customizable and affordable workstation for medium to large Labs. The premium series workstation comes with higher suction capacity and additional features that assists in grossing as well as protects you from inhaling harmful Formalin and other harmful fumes while grossing.



MODEL NO	SIZE (LENGTH)
RB500	5 FEET
RB600	6 FEET
RB700	7 FEET

Specification

Model	RB500
MOC	SS304 (Standard) / SS316 (Optional)
Ventilation	Backdraft
Dimension	1515mm Length x 805mm Width x 2145mm Height
Work Surface	900mm height (Can be customized on request)
Sink	380mm length x 350mm width x 180mm deep sink
Filters	Water washable polyester prefilters along with optional custom designed filter
Electrical	220 - 240 VAC 50/60 Hz, 1 Phase , 110 V supplied on request
Plumbing	1/2" industrial grade flexible hose supplied for water connection 1/2" drain hose (waste disposal)
Waste Disposal	0.55 HP High capacity waste disposal unit.
Exhaust air capacity	1800 CMH
Perforated working area	1150mm length x 460mm width
exhaust duct diameter	200mm diameter

Standard Features



Anti Corrosive Construction



Perforated Worktop



Inbuilt Ventilation



Spacious Sink



Backdraft Extraction



Powerful Lights



Boundry Rinse



Towel Holder




Laser Etched Cutting Board



Instrument Holder (Magnetic)



Custom Filters



Ergonomic Design



Spacious Shelves



Hand Rinse



Bio Decomposer



Transparent Panels



Work floor Levelling



Formalin Dispenser

Optional Features



Magnifying Lamp



Path Imaging System



Safety Shield



Camera Holder



MOC: SS316



Optional Carbon Filter



Pullout Pad



Formalin Monitoring System



Adjustable Shelves



Organizer Bin




Digital Weighing Scale



Ruler (Inch scale)



LED Monitor



CPU | Keyboard Mounting Brackets