
# GROSS IMAGING SYSTEM

## PIXGRAB PRO








**PIXGRAB PRO** brings fast, precise imaging for pathology labs—capture, measure, and annotate gross images with 4K clarity and instant export. Featuring real-time HD video, 35x Zoom, precision autofocus, and user-adjustable workflow controls, PIXGRAB PRO is designed for detail, speed, and documentation compliance. Integrated 3-key pedal switch, Image analysis feature with Measurement and annotation overlay streamlines grossing for every lab.

## Key Features

- 4K Ultra HD Imaging
- JPEG & MP4 Capture
- 1×–35× Digital Zoom
- Auto/Manual Focus
- QR/Barcode Recognition
- LIS Integration (HL7/API/DB)
- Calibrated Measurement
- Image Annotation Suite
- Session-Based Data Management

## Specimen Images







## Technical Specifications

| Parameter | Specification |
|---|---|
| Imaging | 4K live preview, image & video capture. |
| Camera Control | 35× zoom, auto/manual focus & rotation. |
| Measurement | Calibrated distance & area measurement. |
| Annotation | Text, shapes & freehand image markup. |
| Reporting | Undo/Redo with CSV export. |
| Storage | Automatic session-based media storage. |
| Data Management | Custom file & folder naming. |
| Automation | Auto session creation & organization. |
| Connectivity | QR/Barcode scanning & patient identification |
| System Integration | LIS integration via REST API, HL7, DB & CSV. |
| Workflow | Hands-free capture using foot pedal. |
| User Interface | Touchscreen-based intuitive operation. |
| Productivity | Fast, efficient specimen documentation. |

APZEM®

**M: 73388 01069 |** 97890 66648
**E: connect@apzem.com |** sales@apzem.com
**W:**www.grossingstation.com